



University of Wollongong in Dubai

## Lab Experiment: Lenses and Optical Instruments

	1	2	3	4
Family Name:	.....	.....	.....	.....
First Name:	.....	.....	.....	.....
Student Number:	.....	.....	.....	.....

Two Experiments:

- 1) - Lens Equation
- 2) - Galileo telescope

## Experiment 1 Lens equation

### Exercises

- Determine the focal length  $f$  of a lens from the image distance  $b$  and the object distance  $g$ .
- Measure the image distance  $b$ , the image size  $B$ , and the magnification factor  $A$  as functions of the object distance  $g$ .

### Basic principles




Fig. 1: Formation of an image of an object of size  $G$  by a converging lens with a focal length  $f$ .

A converging lens can produce a magnified or a reduced image of an object. The magnification factor  $A$ , which is the ratio of the image size  $B$  to the object size  $G$ , depends on the image distance  $b$  and the object distance  $g$ . The quantities are related by the imaging equation:

$$A = \frac{B}{G} = \frac{b}{g}$$

in which  $G$  is the size of the object to be imaged and  $B$  is that of its image. The object distance and the image distance cannot be varied independently of each other, but are connected by the lens equation:

$$\frac{1}{b} + \frac{1}{g} = \frac{1}{f}$$

### Apparatus required

From basic set for Kröncke optical system (U8477100)

- 1 Optical lamp
- 1 Optical bench
- 4 Optical saddles
- 1 Holding clip
- 1 Screen, white
- 1 Slide with numeral 1
- 1 Lens,  $f = +150$  mm

#### Also required

- 1 Transformer, 12 V, 25 VA (U8475470)
- 1 Ruler

## Experiment set-up




Fig. 2: Slide with numeral 1 in holder, viewed along the direction of the light beam (left) and against the direction of the light beam (right).

- 1) Measure the length of the "1" with the ruler in millimetres and enter the value as the object size  $G$  in the table.
- 2) Fix the slide in a holder with the "1" horizontal and at the centre of the opening in the holder.




Fig. 3: Experiment set-up; side view.

- 3) Position the optical lamp, the holder and the lens with  $f = +150 \text{ mm}$  on the optical bench as shown in Figure 3.
- 4) Connect the optical lamp to the voltage output terminals of the 12-V transformer.

## Experiment procedure

- 1) Set up the screen at the end of the optical bench and slowly move it towards the lens until a sharp image of the "1" is obtained.
- 2) Measure the distance between the slide and the lens and enter the value as the object distance  $g$  in the table (see Fig. 3).
- 3) Measure the distance between the lens and the screen and enter the value as the image distance  $b$  in the table (see Fig. 3).
- 4) Measure the size of the "1" image on the screen with the ruler and enter the value as the image size  $B$  in the table (see Fig. 3).
- 5) By moving the lens, adjust the object distance to  $g = 250 \text{ mm}$  and move the screen until a sharp image of the "1" is obtained.
- 6) Measure the image distance  $b$  and the image size  $B$  and enter the values in the table.

- 7) Repeat the measurements for the other object distances  $g$  listed in the table.
- 8) Calculate the ratios  $B/G$  and  $b/g$  and enter the results in the table.
- 9) From the measured values of the object distance  $g$  and the image distance  $b$ , use the lens equation to calculate the focal length  $f$  and enter the results in the table.
- 10) From the calculated focal length values, determine the average value of  $f$ , and set the object distance  $g$  to  $2f$ , twice the focal length.
- 11) For the latter special case, measure the image distance  $b$  and the image size  $B$  and enter the values in the table.

Data table:

$G = 23 \text{ mm}$

$g / \text{mm}$	$b / \text{mm}$	$B / \text{mm}$	$\frac{b}{g}$	$\frac{B}{G}$	$f / \text{mm}$
200					
250					
300					
450					
600					
Average of focal length values:					
$2f =$					

### Review and evaluation of results

- a) Describe the appearance of the image of the "1" on the screen.
- b) How does the image size depend on the image distance?
- c) How does the image size depend on the object distance?
- d) How are the ratios  $B/G$  and  $b/g$  related to each other?
- e) What is the focal length of the lens?
- f) For which object distances is the image reduced and for which is it magnified?
- g) When are the image and the object of equal size?

## Experiment 2 Galileo telescope

### Exercises

- Construct a model of a Galileo telescope.
- Determine the angular magnification and the required distance between the lenses as functions of their focal lengths.

### Basic principles




Fig. 1: How a Galileo telescope magnifies the angle of vision  $\epsilon_0$  for a distant object.

When a far distant object is viewed through a Galileo telescope it appears larger, because the telescope magnifies the angle of vision that the object subtends at the eye of the observer. The angular magnification is achieved by using a converging lens and a diverging lens separated by a distance  $d$  that is exactly equal to the sum of their focal lengths, taking into account the fact that the focal length of any diverging lens is treated as negative:

$$d = f_1 + f_2$$

The converging lens, the objective, gives rise to a real inverted image of the distant object in its focal plane, with size  $B$  given by:

$$B = f_1 \cdot \epsilon_0$$

This image is then viewed with the eye focused on infinity through the diverging lens, the ocular. The resulting increased angle of vision is given by:

$$\epsilon = \frac{B}{|f_2|} = \frac{f_1}{|f_2|} \cdot \epsilon_0$$

Thus, the original angle of vision  $\epsilon_0$  subtended at the naked eye by the object has been magnified by the factor:

$$A = \frac{f_1}{|f_2|}$$

### Apparatus required

From basic set for Kröncke optical system (U8477100)

- 1 Optical bench
- 2 Optical saddles
- 1 Lens,  $f = +150$  mm
- 1 Lens,  $f = +300$  mm
- 1 Lens,  $f = -100$  mm

## Experiment set-up




Fig. 2: Experiment set-up; side view.

Set up the  $-100\text{ mm}$  lens (the ocular) and the  $+300\text{ mm}$  lens (the objective) on the optical bench as shown in Figure 2.

## Experiment procedure

Use the telescope model to view far distant objects such as houses or trees, or – all the better for estimating the angular magnification – a distant wall-mounted ruler divided into two differently coloured parts.

- 1) Look through both lenses at the distant object, and obtain a sharp image by moving the objective lens.
- 2) Measure the distance  $d$  between the ocular and the objective and enter the value in the table.
- 3) Calculate the sum of the focal lengths,  $f_1 + f_2$ , and compare it with the distance  $d$  between the lenses.
- 4) Look at the telescope image with one eye and look past the lenses at the object itself with the other eye.
- 5) Compare the directly viewed object with its magnified image, estimate the angular magnification  $A$ , and enter the value in the table.
- 6) Calculate the ratio  $\frac{f_1}{|f_2|}$  and compare it with the estimated angular magnification  $A$ .
- 7) Repeat the experiment with the  $+150\text{ mm}$  lens instead of the  $+300\text{ mm}$  lens.

Data table:

$f_1 / \text{mm}$	$f_2 / \text{mm}$	$d / \text{mm}$	$f_1 + f_2 / \text{mm}$	$A$	$\frac{f_1}{f_2}$
150	-100				
300	-100				

### Review and evaluation of results

- a) Describe the size and appearance of the image seen through the ocular.
- b) What is the distance between the lenses that gives a sharp image?
- c) What is the relationship between the angular magnification  $A$  and the ratio  $\frac{f_1}{|f_2|}$ ?
- d) Which of the objective focal lengths gives the greatest angular magnification?
- e) What is the relationship between the angular magnification and the length of the telescope?

